



**PrimaLaport Optical, PrimaLaport Bladed, PrimaLaport Sleeve,
PrimaLaport Hasson, PrimaLaport Mini port bladeless,
PrimaLaport U-Flex**

Description

Sterile Trocars for Endoscopic surgery, PrimaLaport, is a sterile single use device and is intended to establish a path of entry for endoscopic instruments for use during general, abdominal, gynaecological and thoracic minimally invasive procedures or to gain access through tissue planes and/or potential spaces for endoscopic instruments.

■ Rotating reducer cap; It has a reducer cap, which can be rotated 360°, for a comfortable performance of surgical procedures.

■ One touch valve control; A valve, which is manipulated with a button, allows sealing of the system. Its activation is useful for intra- abdominal gauze or specimen retrieval. It can be activated with a single hand. A stopcock for gas insufflations is included.

■ No fear for reducer cap loss; The reducer cap is connected to the main body of trocar through a hinge, thus risk of reducer cap loss or contamination is eliminated.

■ Simple blade activation system; The obturator has a safety button which activates a retractile internal cutting knife. It is only activated during trocar insertion maneuver.

Product Specifications

Please refer to 'description' on the Inner Box.

Intended Purpose

It is a needle and/or sleeve that is used for the centesis or tresis of tissue in order to insert medical device into human body during endoscope surgery and is disposable.

Trocar Laparoscópico estéril de uso único Instrução para Uso

Indications

- Used in laparoscopic surgery
- Used in thoracoscopy surgery

Cautions during use



Warnings

- 1) User shall consider the position of operated area when selecting device.
- 2) User shall consider the size of device in order to eliminate the occurrence of complication.
- 3) If you have a metal allergy, please consult your doctor.
- 4) If the device is reused, it may cause to patients infection.
- 5) Do not reuse, reprocess or re-sterilize this device.
- 6) Discard once released product though disused.
- 7) Do not use and discard when sterile packaging is damaged. If the packing of product is broken, it may cause to patients infection.
- 8) Do not use after expiration date.

Cautions

- 1) Do not use with damaged or defective device.
- 2) This product has to be used by specialist who has surgical operation experience of laparoscope and so on.
- 3) Remove the product from sterile package just before using and use it. If time is elapsed on opened condition, product may be contaminated and patient may be infected.
- 4) Avoid the direct rays of the sun, store at temperature range (2~30°C).
- 5) Be careful not to damage sterile package on storage and transfer.

Instruction For use

Optical trocars (PrimaLaport Optical)

1. Take out PrimaLaport from the blister package.
2. Insert trocar needle in head of trocar sleeve.
3. Insert endoscopic camera in hole of trocar needle.
4. Puncture beforehand the abdominal cavity with surgical blade on surgical region.
5. After inserted blade of trocar needle in the abdominal cavity, trocar sleeve push carefully when tip of trocar sleeve passes the abdominal cavity.
6. Remove trocar needle, regulate trocar sleeve's depth according to user's the convenience and use.
7. If discharge smoke and gas, press trocar button.
8. Discard this product after using because this is disposal.


Bladeless trocars (PrimaLaport Mini port bladeless, Prima Laport U-Flex, PrimaLaport Hasson)

1. Take out PrimaLaport from the blister package.
2. Insert trocar needle in head of trocar sleeve.
3. Puncture beforehand the abdominal cavity with surgical blade on surgical region.
4. After inserted blade of trocar needle in the abdominal cavity, trocar sleeve push carefully when tip of trocar sleeve passes the abdominal cavity.
5. If using the PrimaLaport Hasson, bind skin and wings of Hasson using the tread.
6. Remove trocar needle, regulate trocar sleeve's depth according to user's the convenience and use.
7. If discharge smoke and gas, press trocar button.
8. Discard this product after using because this is disposal.

Bladed trocars (PrimaLaport Bladed)

1. Take out Laport from the blister package.
2. Push button at handle of trocar needle.
3. Insert trocar needle in head of trocar sleeve.
4. Puncture beforehand the abdominal cavity with surgical blade on surgical region.
5. After inserted blade of trocar needle in the abdominal cavity, trocar sleeve push carefully when tip of trocar sleeve passes the abdominal cavity.
6. Remove trocar needle, regulate trocar sleeve's depth according to user's the convenience and use.
7. If discharge smoke and gas, press trocar button.
8. Discard this product after using because this is disposal.

Use, Storage and Transport Conditions

	<ul style="list-style-type: none">- Please store this device under a low temperature (2~ 30°C).- Avoid the direct sunlight.
---	--


Compatibility Device

The device can be available in all products with laparoscopic trocar or singleport on the market.

However, the size of the connected with the luer-lock insufflation valve or connectors shall be Ø4.25.

Shelf life : 5Years

Symbols Used for Labeling

	Manufacturer		Model name
	Consult instructions for use		Date of Manufacturing
	Lot Number		Caution
	Temperature limitation		Keep away from sunlight
	Use by date		Sterilized using ethylene oxide
	Do not Use if package is damage		Do not reuse
	Handle with care		Keep dry
	Fragile		Latex free

Distributed by:

LAWMED

Unit 1, Walton Lodge, Bridge Street
Walton On Thames, KT12 1BT
United Kingdom
ORDER LINE +44 (0)1932 260 838

Manufactured by:

SEJONG MEDICAL Co.,Ltd.
11, Shinchon 2-ro, Paju-si Gyeonggi-do,
10880, South Korea
TEL: +82-31-945-8909
EMAIL: info@sejongmedical.com



Obelis S.A
Bd Général Wahis, 53
B-1030 Brussels, Belgium