

Lapsingle Search Port Wound Retractor



Sterile Single use Laparoscopic port and Wound retractor
Instructions For Use

Lapsingle / Lapsingle Plus / Lapsingle Vision /
Search Port / Wound Retractor

Description

This is single or multiple access port used in single incision laparoscopic surgery which is combined with wound retractor. It spreads the incision site into an annular shape in order to use multiple instruments while blocking gas emissions

Product Specifications

Please refer to 'description' on the Inner Box.

Intended Purpose

This is intended to establish a path of entry for laparoscopic instruments for use during minimally invasive abdominal laparoscopic surgery.

Cautions during use



Warnings

- 1) Do not use the device to the patient who is not applicable.
- 2) This product is sterilized with EO gas, intended for single use and should not be reprocessed or re-sterilized.
- 3) Re-sterilization or reuse may compromise the integrity of the device and may cause hazardous injury or infection to the patients.



Cautions

- 1) The product must be maintained in sterilized condition and do not use the product if the term of validity is expired.
- 2) Do not use if the sterile package is open or damaged.
- 3) Do not reuse, reprocess or re-sterilize this device.
- 4) Do not disassemble or modify the device without permission of the manufacturer.
- 5) Do not use besides from intended use.
- 6) Do not use if there is any damage on the package or product.

Precaution For use

This is a sterile product, check the package and the expiration date prior to use. Device should be performed only by the trained physicians. Read the instruction for use carefully prior to use.

Instruction For use

Lapsingle, Lapsingle Plus, Lapsingle Vision

1. Remove the product from the package.
2. Make an incision about ±20mm at the umbilical.
3. Insert the lower ring of wound retractor in the incision site.
4. Expand the incision site by rolling the wound retractor and fix it to the incision site by tightening the upper ring.
5. Assemble Lapsingle with the upper ring of wound retractor.
6. Insufflate the gas through the luer-lock insufflation valve.
7. Use the right size port for surgical procedure.
8. When the surgical procedure is done, disassemble the Lapsingle with the upper ring of wound retractor.
9. STRING TYPE : Pull the string of wound retractor after the surgical procedure for removal.
Non-STRING TYPE : Pull the wound retractor after the surgical procedure for removal.

Lapsingle Tube, Lapsingle Plus Tube, Lapsingle Vision Tube

1. Remove the product from the package.
2. Make an incision about ±20mm at the umbilical.
3. Insert the lower ring of wound retractor in the incision site.
4. Expand the incision site by rolling the wound retractor and fix it to the incision site by tightening the upper ring.
5. Assemble Lapsingle with the upper ring of wound retractor.
6. Insufflate the gas through the luer-lock insufflation valve.
7. Use the right size port for surgical procedure.
8. When the surgical procedure is done, disassemble the Lapsingle with the upper ring of wound retractor.
9. STRING TYPE : Pull the string of wound retractor after the surgical procedure for removal.
Non-STRING TYPE : Pull the wound retractor after the surgical procedure for removal.

Search Port


1. Remove the product from the package.
2. Make an incision about ±20mm at the umbilical.
3. Insert the wound retractor to the incision site.
4. Roll the spreader to spread the incision site.
5. Place the cover of Search Port in incision site and assemble it.
6. Insufflate the gas into the abdominal and check through the Search port whether the surgical site is well-finished.
7. Disassemble the Search port with the upper ring of wound retractor.
8. STRING TYPE : Pull the string of wound retractor after the surgical procedure for removal.
Non-STRING TYPE : Pull the wound retractor after the surgical procedure for removal.

Wound Retractor

1. Remove the product from the package.
2. Make an incision about ±20mm at the umbilical.
3. Insert the lower ring of wound retractor in the incision site.

4. Expand the incision site by rolling the wound retractor and fix it to the incision site by tightening the upper ring.
5. STRING TYPE : Pull the string of wound retractor after the surgical procedure for removal.
Non-STRING TYPE : Pull the wound retractor after the surgical procedure for removal.

Use, Storage and Transport Conditions



Biocompatibility

Sterile Laparoscopic single port and wound retractor is made by the materials which are confirmed to be biocompatible in compliance with ISO 10993 series standards.

Shelf life : 3 Years

Symbols Used for Labeling

	Manufacturer		Model name
	Consult instructions for use		Date of Manufacturing
	Lot Number		Caution
	Temperature limitation		Keep away from sunlight
	Use by date		Sterilized using ethylene oxide
	Do not Use if package is damage		Do not reuse
	Handle with care		Keep dry
	Fragile		Latex free
	Korea Good Manufacturing Practice		CE Mark
	Warning		EU Authorized Representative



Manufactured by Sejong Medical Co.,Ltd.
11, sinchon 2-ro, Paju-si, Gyeonggi-do, 10880, South Korea
Tel: +82-31-945-8909 Fax: +82-31-945-8190
www.sejongmedical.com

OBELIS S.A
Boulevard General Wahis 53 1030 Brussels, BELGIUM
SJ-PG-21-IFU02 (Rev.2.1_2018-07)