

LapShieldTM

Wound Retractor

Elastic Abdominal Retractor for Minilaparotomy

- Maximize exposure, minimize incision size
- Reduce the risk of wound infections
- Offer ultimate versatility



360° Atraumatic Retraction

- Distributes force evenly for optimal field of view
- Allows for maximum exposure with a minimum incision size
- Soft materials and evenly dispersed tension promote patient post-operative comfort

360° Protection


- Significantly decreases wound infection
- Allows visualization of wound margins
- Maintains moisture at the incision site
- Helps protect wound from contamination while helping to keep wound edges moist

Ultimate Versatility

- Unique design facilitates rapid and effortless set-up
- Multiple sizes accommodate a wide range of surgeries

Easily Placed and Removed

- Inserted in seconds
- Flexible ring slips gently into incision
- Removed with one finger



Product Code	Incision Length (mm)	Quantity
* WR-30US	10~20	6 /box
WR-60ES	20~40	6 /box
WR-80S	30~60	6 /box
WR-110M	50~90	6 /box
WR-160L	90~140	6 /box
WR-190EL	110~170	6 /box

* Availability to be advised

Designed for soft tissue retraction, LAGIS LapShield™ provides unparalleled exposure without the trauma associated with prolonged point retraction.