



MID-TUBE[®]

OROGASTRIC CALIBRATION TUBE



36 Fr



38 Fr



MID[®]
INNOVATION FOR SURGERY

MID-TUBE[®] FEATURES:

MID-TUBE[®] is a silicone tube with a diameter of 36 Fr or 38 Fr.

MID-TUBE[®] can also be connected to a suction system for easy removal of air or fluid in the patient's stomach either before or during the procedure (leak test).

6 graduations (from 35 to 50 mm).

EXAMPLES OF PROCEDURE IN OBESITY SURGERY:

Gastric by-pass
Sleeve gastrectomy

...



TECHNICAL DATA:

Material used:	Medical Silicone
Tube length:	800 mm
Tube external diameter:	12 mm (36Fr) 12.8 mm (38Fr)

PACKAGING & PRODUCT CODE:

1 unit packaged in a sealed pouch per box.

Shelf box dimensions: 24.6 x 17.5 x 3.4 cm.

REF **MID130**

REF **MID136**

WARNINGS:



The MID-TUBE[®] must only be handled by practitioners with experience of digestive surgery.



The particular instructions for use are shown in the notice.



Single Use / Do not resterilize.



Contents are sterile if package is unopened and undamaged.

CE 0459 - Classe: IIa



Sterile (Ethylene Oxide).



Store in a cool, dark, dry place.

MID 
INNOVATION FOR SURGERY

VOTRE CONTACT

9 CHEMIN DU JUBIN
69570 DARDILLY
FRANCE

Tel.: +33 (0)4 78 17 48 04

Fax.: +33 (0)4 72 82 91 23

mid@mid-med.com



www.mid-med.com

This documentation is only intended for healthcare professionals. This document is the technical file of MID13X MID-TUBE[®].