



PRIMA·bag  
Laparoscopic  
retrieval pouch



ADVANCING SURGICAL POSSIBILITIES

# PRIMA<sup>•</sup>bag

The Prima bag is designed with a reinforced seam and flexible pouch to ensure the safe and easy removal of tissue within the abdominal cavity during a laparoscopic procedure

The Prima bag has an automatic mechanism and is available in a wide range of sizes to accommodate all tissue types. The range extends from a small Y-type bag to an extra large U-type bag

All Prima bags are compatible with a 10 mm trocar except for the XX-large which enters the abdominal cavity through a 12 mm trocar, removing the need to open an additional 15 mm trocar

Category	Length	Width	Volume	Trocar size	Units per box	Product codes
Y-TYPE	140 mm	70 mm	75 ml	10 mm	5	B21-1032
V-TYPE	140 mm	100 mm	225 ml	10 mm	5	B21-1033
U-TYPE, Medium	125 mm	125 mm	410 ml	10 mm	5	B21-1003
U-TYPE, Large	160 mm	160 mm	920 ml	10 mm	5	B21-1004
U-TYPE, X-Large	200 mm	200 mm	1540 ml	10 mm	5	B21-1005
U-TYPE, XX-Large	300 mm	250 mm	4800 ml	12 mm	5	B11-1208

## ADVANCING SURGICAL POSSIBILITIES

LAWMED Ltd.  
Unit 1, Walton Lodge, Bridge Street  
Walton On Thames, KT12 1BT  
United Kingdom

E: [enquiries@lawmed.co.uk](mailto:enquiries@lawmed.co.uk)  
T: +44 (0)1932 260 838  
F: +44 (0)1932 260 849  
[www.lawmed.co.uk](http://www.lawmed.co.uk)