




Advancing
Surgical
Possibilities

Introducing the **iScalpel™**

The Prima iScalpel enables precise open and laparoscopic dissection using ultrasonic energy to cut and coagulate tissues and vessels up to 5mm in diameter

Item description	Instrument length	Units per box	Product codes
5mm ultrasonic multi-purpose shears	36cm	3	USE36
5mm ultrasonic multi-purpose shears	23cm	3	USE23
5mm ultrasonic multi-purpose shears	14cm	3	USE14
Reusable hand piece	N/A	1	HP25



For more information please contact us at:

E: enquiries@lawmed.co.uk

T: +44 (0)1932 260 838

www.lawmed.co.uk

LM37 SEPT 2018