



PRIMA•laparoscopic

Item Name	Item Description	Units Per Box	Product Codes (330 mm Length)	Product Codes (450 mm Length)
Prima babcock	Babcock with ratchet	6	PMB-1	-
Prima clinch	Clinch style grasper with ratchet	6	PMC-1	PMC-145
Prima dissect	Maryland dissector	6	PMD-1	PMD-145
Prima grasp	Maryland dissector	6	PMG-1	PMG-145
Prima johan (22cm jaw)	Atraumatic grasper with ratchet	6	PMJ-1	PMJ-145
Prima johan 2 (22mm jaw)	Johan style grasper with ratchet	6	PMJ-2	-
Prima metz (insulated to tip)	Metzenbaum scissor	6	PMM-1	PMM-145

All devices include monopolar diathermy

The Prima laparoscopic range is designed with versatility in mind. 7 different tip configurations ensure the correct instrument is available for each procedure regardless of its complexity.

Each instrument tip is easily identified by the unique colour on its rotatable dial, simplifying instrument selection during a procedure.

Prima laparoscopic 5 mm instruments are available in 2 lengths, 330 mm and 450 mm.

