


## KOL Cecil – Dedicated for TaTME Cecil Approach



## Intraluminal Stapler for Single Use

Product Code	Anvil Diameter	Colour Code	Staple Quantity	Blade Diameter	Staple Height	Closed Staple Height	Housing Length / Volume
CSC29-KOL	29 mm	Blue	24	20.5 mm	5.0 mm	1.0 ~ 2.5 mm	4 cm / 12.6 cm <sup>3</sup>
CSC33-KOL	33 mm	Green	32	24.8 mm	5.0 mm	1.0 ~ 2.5 mm	4 cm / 17.6 cm <sup>3</sup>

Transanal Accessories



Reference:  
 Ma et al. BMC Cancer (2016) 16:380 Doi: 10.1097/sla.0000000000003757  
 Lacy et al. BMC Cancer (2020) 20:677  
 Surgical Endoscopy 33,2459–2467(2019)



ADVANCING  
 SURGICAL  
 POSSIBILITIES  
 TOGETHER

Distributed by  
 Lawmed Ltd.  
 Unit 1, Walton Lodge, Bridge Street  
 Walton On Thames, KT12 1BT  
 United Kingdom  
 E: enquiries@lawmed.co.uk  
 T: +44 (0)1932 260 838  
 www.lawmed.co.uk

**KOL Cecil**  
 Dedicated for TaTME Cecil Approach

# KOL Cecil

## Dedicated for TaTME Cecil Approach



### Bigger Anvil Housing


Easier purse string knotting  
Better anastomosis formation  
by allowing more tissue  
inside the recess

### Long Spike


Optimizing for TaTME

### Suture Hole

Easier transanal  
manipulation  
for TaTME



### Double Space Specimen Housing



### Four Traction Holes

Allow for the pulling of  
the staple line and tissue  
into the housing