


Intraluminal stapler for single use



Intraluminal stapler for single use



Product code	Colour code	Staple quantity	Anvil diameter	Blade diameter
ECSC21	Orange	16	21mm	13.0mm
ECSC25	White	20	25mm	17.0mm
ECSC29	Blue	24	29mm	20.5mm
ECSC33	Green	28	33mm	24.8mm
CSC21-A	Orange	16	21mm	13.0mm
CSC25-A	White	20	25mm	17.0mm
CSC29-A	Blue	24	29mm	20.5mm
CSC33-A	Green	28	33mm	24.8mm
CSC21	Orange	16	21mm	13.0mm
CSC25	White	20	25mm	17.0mm
CSC29	Blue	24	29mm	20.5mm
CSC33	Green	28	33mm	24.8mm



Manufactured by:

TOUCHSTONE

Touchstone International Medical Science Co., Ltd. 21A Science Plaza, International Science Park, No.1355 Jinjihu Avenue, Suzhou, 215021, P.R.China
Tel: +86 512 6299 1933; 6299 1939

Protected by one or more of the following patents:

200710191747.X 200710190987.8 200710191026.9 200710192011.4
200710190965.1 200520074202.7 200520072023.X 200620076828.6
200610041509.6 200820031086.4 200510038780.X 200520068606.5
200520070443.4 200520074201.2 200520074204.6 200520139870.3
200520139728.9 200520139783.8 200520139784.2 200520139726.X
200620073826.1 200620125458.0 200520139868.6 200620125462.7
200620078175.5 200710191957.9



ADVANCING SURGICAL POSSIBILITIES

LAWMED Ltd.

Unit 1, Walton Lodge, Bridge Street
Walton On Thames, KT12 1BT
United Kingdom

E: enquiries@lawmed.co.uk
T: +44 (0)1932 260 838
www.lawmed.co.uk

ADVANCING SURGICAL POSSIBILITIES

Endoscopic circular stapler curved




Unique rotating blade system

Clinical benefits

RBS™ : Advanced technology beyond the conventional cutting to reduce leakage and bleeding.

"From the 100 bariatric patients, we have had zero leaks and zero bleeding"
 - By Dr. Rudolf Steffen, Switzerland Reference: *"Time for a change in surgical stapling - introducing the Touchstone Circular Staple gun"* by Dr. Rudolf Steffen, Center of Excellence Hospital for Bariatric Surgery Beau-Site Clinic, Berne, Switzerland.

Bariatric News, Issue 4, March, 2010



ASL™ : Auto Safety Lock System to prevent 2nd misfiring and convey the audible feedback to the surgeon.

Reference: Used in the laparoscopic By-pass by Prof. Karl Miller, head of Surgical Department, Hallein Clinic, Austria.



Bendos™ Hole: Facilitate insertion of gastric feeding tube to avoid damaging the anastomosis.

