


Advancing  
surgical  
possibilities



LINEAR STAPLER AND RELOADS

# Linear Stapler and Reloads

*Reloadable 11 times  
for a total of 12 firings*




2-Step Controlled Closing

Release Button

Clear Firing Indicator

To Close Jaws



To Fire



To Open



## Unique Features

- 2-Step closure - allows surgeon to adjust the tissue before closing the jaws completely
- Single hand operation - freeing other hand for any assistance
- Locking handle - provides clear firing indicator
- Choice of manual or automatic pin placement
- Cartridge safety plate



## Linear Stapler and Reloads

Product Code	Staple Line Length	Rows of Staple	Staple Quantity	Staple Height	Colour Code
LS3035	30	2	13	3.5	Blue
LSC3035	30	2	13	3.5	Blue
LS3048	30	2	13	4.8	Green
LSC3048	30	2	13	4.8	Green
LS30V3	30	2	23	2.5	White
LSC30V3	30	2	23	2.5	White
LS6035	60	2	25	3.5	Blue
LSC6035	60	2	25	3.5	Blue
LS6048	60	2	25	4.8	Green
LSC6048	60	2	25	4.8	Green
LS9035	90	2	33	3.5	Blue
LSC9035	90	2	33	3.5	Blue
LS9048	90	2	33	4.8	Green
LSC9048	90	2	33	4.8	Green

**Advancing surgical possibilities**

LAWMED Ltd.  
Unit 1, Walton Lodge, Bridge Street  
Walton On Thames, KT12 1BT  
United Kingdom

E: [enquiries@lawmed.co.uk](mailto:enquiries@lawmed.co.uk)  
T: +44 (0)1932 260 838  
F: +44 (0)1932 260 849  
[www.lawmed.co.uk](http://www.lawmed.co.uk)

Manufactured by  
**TOUCHSTONE™**

Distributed by  
LAWMED Ltd.